

coronary guidewire
Conductor
Precision at Your Fingertips




concertmedical

Concert Medical, LLC • 77 Accord Park Drive • Norwell, MA 02061
Ph: 781-871-7882 • Fax: 781-871-6657 • www.concertmedical.com

The **Conductor** Guidewires are designed to provide
**Consistent Performance and Support, Distinct 1:1 Torque,
 Tip Shapeability and Tip Shape Retention.**

CONDUCTOR GUIDEWIRE FAMILY:
 0.014" x 190 cm / 300 cm



Conductor Very Floppy




Consistent Performance During the Most Demanding Cases
 Most Atraumatic – suitable for sensitive occlusions
 Tip flexibility: 0.90 g
 Support: 40 g/cm (10 cm out, at 15° bend)
 Order Number: 511950-001

Conductor Floppy



Workhorse Guidewire with Precise Torque Control and Easy Tip Shapeability
 Great durability and tip shape retention
 Tip flexibility: 1.3 g
 Support: 49 g/cm (10 cm out, at 15° bend)
 Order Number: 511950-002

Conductor Medium



3 Gram Tip Design with Moderate Support
 For complex cases and tight lesions
 Tip flexibility: 3g
 Support: 58g/cm (10 cm out, at 15° bend)
 Order Number: 511950-003

Conductor Extra Support



Designed for Optimal Performance and Deliverability
 Recommended for direct stenting and complex anatomies
 Tip flexibility: 1.3 g
 Support: 112 g/cm (10 cm out, at 15° bend)
 Order Number: 511950-006

Also available: **Concert EX** - Extension Guidewire
 0.014" x 152 cm. Order number: 511920-001

signature